

Service Data Tool 2 for OpenBus™ Control Boxes and Batteries

User Manual



Contents

Introduction	4
Usage	4
Before getting started with reading out service data on a PC	5
OneConnect	6
SDT2 for CO control boxes	7
System overview for CO control boxes	7
Reading out service data on a PC	8
User interface	9
Summary view	10
Actuator view	12
How to conclude on the service information	13
Controls view	14
Report view	15
Advanced settings view	17
Troubleshooting view	18
PJ2 Junction box view	19
Battery view.....	20
BA22 Battery Service view.....	21
BA16/19 Battery Service view.....	22
How to change the application picture, company logo and information PDF-file	23
SDT2 for CO-Link systems	24
System overview for CO-Link systems	24
CO-Link system explanation.....	25
SDT2 for LIFT50 systems	26
System overview for LIFT50.....	26
Reading out service data on a PC.....	27
LIFT50 Service Data Tool sections	28
Summary view	29
Actuator view	31
COL50 Controls view.....	33
COL50 Report view.....	34
COL50 Advanced settings.....	37
BAL50 Battery Service View	40
How to change the application picture, company logo and information PDF-file	41

SDT2 for JUMBO Care systems	42
System overview for JUMBO Care	42
Service intervals	44
Read out service data on a JUMBO Care with display	45
Reading out service data on a PC.....	46
Actuator view	47
How to conclude on the service information	48
Control box view	48
Controls view	49
Report view	50
Advanced settings view	51
Resetting of service interval after service has been carried out.....	52
Info site	52

Introduction

Service Data Tool ensures efficient maintenance of healthcare applications

Service Data Tool enables easy read-out of service data from OpenBus™ control boxes and batteries.

Usage

- Compatible with OpenBus control boxes, CB6 MK2, BA16, BA19, BA2, CBJC, BAL50 and more
- Compatible with Windows 10 and 11



Batteries supported in SDT	PCP 1.0	PCP 2.0
BA16		X
BA19		X
BA22		X
BAL50		

The version number on the product label will tell the PCP version

If the version number is 1.xx or less, the battery is PCP 1.0

If the version number is 2.xx or more, the battery is PCP 2.0

Before getting started with reading out service data on a PC

Equipment needed to read out service data

When reading out service information on a PC you need:

- Service Data Tool 2 version 2.7.9 or newer version installed on the PC.
- The software will be provided from your local LINAK sales representative.
- An OpenBus™ programming and data read out box (LINAK item number IB300001)
- A service readout cable (LINAK item number 0964198)
- Modular junction box for connection of Service Data Tool (LINAK item number MJB000(3/4/5)000-1023)
- Cable for connection of the Modular junction box to COXX (LINAK item number 0964461-XXXX-A)
- One USB A-B cable

How to connect the equipment and get started

1. First, make sure that Service Data Tool drivers and the Service Data Tool 2 software (version 2.7.9 or higher) is installed on your PC. If Service Data Tool drivers and the Service Data Tool software are not installed, please contact your local LINAK subsidiary for the install files. Connect the equipment as shown in the following system overviews.
2. Wake the COXX (by pressing a button on the hand control or by pressing the foot switch).
3. Enter the Start menu on your PC and open the Service Data Tool (SDT2).
4. When connection is successful, the red indicator ("Supply 8V Missing") in the current PC window changes to green.



OneConnect



For a simple, wireless alternative to Service Data Tool that can be used with Android and iOS devices, please see OneConnect by LINAK.

OneConnect is a powerful app designed to better connect you to your end users and improve troubleshooting like never before.

- View system status on the fly
- Access system data via Bluetooth®
- Remotely access system data from anywhere in the world
- Customise your help content
- Brand to your company's aesthetics

To learn more about these features and how to get started with OneConnect for your products, visit the [OneConnect landing page](#) and contact us today.

Compatible LINAK Bluetooth control boxes are: COL50, CO71, CO61 (select versions), CA63, CAL40+, and CO53 (select versions).




SDT2 for CO control boxes

System overview for CO control boxes

The LINAK IB3 tool (the SDT “box”) allows for service data read-out from most LINAK systems with a modular (also called RJ50) connection. Below is a diagram of how this IB3 tool is often connected to a standard LINAK system via a Modular Junction Box (MJB).

While it is also possible to connect the yellow connector directly to the CO control box, it is advisable to do so through an MJB so other components can still be connected and ready to send commands, both to wake up the system for the SDT2 program, but also to run the system for troubleshooting purposes.



Information


Please remember to order cables, the modular junction box and the IB300001 box for data read-out from COXX to PC.

Article number	Product
IB300001	Box for data read out from COXX primary or COXX secondary to PC
MJB000 (3/4/5) 000-1023	Modular Junction box for connection of Service Data Tool and foot switch, hand control or other OpenBus™ accessories
0964198	Cable for usage of Service Data Tool (from MJB to IB300001 box)
0964461-xxxx-A	Cable for connection of MJB to COXX: Straight, RAL 7035, Modular Jack, RAL7035, Modular jack, length xxx mm
SML912497-A	Splitter mains cable for COXX primary and COXX secondary




Reading out service data on a PC

Ensure that the COXX or battery view is initiated by pressing the menu shown.
 Please contact your local LINAK Supplier for support if this is not the case.




Indicates that the connection is working 100%

Indicates something is missing or has an error. A small text right to the LED informs about the problem.



Connection to IB30 box OK

Control Box connected = 8V missing











Com Port Connection SDT IB30 version


OpenBus Ok/Fail




User interface







The Service Data Tool for COXX is developed with focus on user friendly design which is easy to use. Information is divided into sections as below:

	Symbol	Explanation
A		For quick and easy overview of actuators info and fatal error status. Please notice: A new actuator connected is only visible in the Service Data Tool view if it has been running in outward or inward direction – A*S>0.
B		For detailed information about the actuators' statistical service data. For refreshing data from the control unit connected. For update of information when replacing the actuator.
C		For detailed information about hand control signals and codes.
D		For service reporting, production number, software number, item number of the application. Saving complete and relevant information for actuators data and service.
E		Intended for trained and authorised service technicians only. For update of actuator info if the control is replaced to maintain service data.
F		For detailed information about errors on the actuator, control box or the controls. Easy download of log file to send for further technical support
G		Please notice! This section is only visible, if the control box software specifies actuators used on Port Junction box PJ2. For detailed information about Port Junction Box
H		Please notice! This section is only visible when COXX is connected to LINAK Battery lithium Ion battery. For information about battery life, charging state and use.
J	Status bar	Status is either IDLE or Reading Extended OpenBus device 8 register 20 It means that SDT is reading data – please wait.





Sections are updated every second to keep data updated. Please re-start SDT if the control box has powered down or disconnected. For further help and recommendations on each section, please press the help icon 

Summary view




Symbol	Explanation
<p>A</p> 	This symbol indicates that COXX is connected to mains power.
	<p>This symbol indicates that COXX is operating on a Lithium Ion battery.</p> <p>Battery Status is indicated by % and colour:</p> <p>Green: Battery fully charged (approx. 100-40% capacity remaining)</p> <p>Yellow: Charging recommended (approx. 40-21% capacity remaining)</p> <p>Orange: Low critical battery level. (Depending on the COXX SW there will be audio signal when the hand control is activated. Limited the actuator function).</p>
	<p>This symbol indicates that COXX is operating on a BA19 battery, Lead Acid.</p> <p>Regularly charging or charging before use is recommended as battery status indication is not available.</p> <p>Green: Battery will remain green. Please notice! There is no low battery warning.</p>
<p><input type="radio"/> Duty cycle guard</p>	If the LED is orange, it indicates that the COXX has stopped operating due to the duty cycle guard.
<p><input type="radio"/> Current fatal error</p>	If the current fatal error LED is orange then you have fatal error and can use the troubleshooting section and button to for help.
<p>D</p> <p>ACT 6 <input checked="" type="checkbox"/></p>	Selected or de-selected the actuators you want to monitor
<p>ACT 6 <input checked="" type="checkbox"/></p> 	The actuator position is shown with numbers.
	<p>The green arrow indicates that the actuator is working and shows which direction, which is also written.</p> <p>Status can be; Direction In, Direction Out, EOS in, EOS out or IDLE</p>
<p>ACT 6 <input checked="" type="checkbox"/></p> 	If the actuator has lost its position, POS LOST indicator will occur in the status bar






	Symbol	Explanation
E		Help manual of what the different functions means.
F		The customer can add their own information as PDF-file.
G		A picture of the control box
H		The customer can change the application picture.
I	Company logo	The customer can add their own LOGO.






Actuator view



Example given with a 6-actuator system

Symbol	Explanation
<p>A</p> <p>Cutoff Limit</p> <p>IN OUT</p> <p>[A]</p> <p>3 3</p>	<p>To ensure that the application will stop if the current draw exceeds the preset limits.</p>
<p>Work</p> <p>A * S</p> <p>21</p>	<p>Total work on the actuator (A*S): Work indicator for the actuators measures via ampere usage *seconds in use. The work indicator gives a very good indication of how much the actuator is worn.</p> <p>LA20: 10.000 cycles in life test equals: 1.300.000 A*S</p> <p>LA23: 10.000 cycles in life test equals: 1.400.000 A*S</p> <p>LA27: 10.000 cycles in life test equals: 3.700.000 A*S</p> <p>LA31: 10.000 cycles in life test equals: 3.000.000 A*S</p> <p>LA34: 10.000 cycles in life test equals: 4.300.000 A*S</p> <p>LA43: 10.000 cycles in life test equals: 2.500.000 A*S</p> <p>LA40: 10.000 cycles in life test equals: 3.800.000 A*S (LA40 6000 N)</p> <p>LA40: 10.000 cycles in life test equals: 4.400.000 A*S (LA40 PL 8000 N)</p> <p>LA40: 10.000 cycles in life test equals: 2.900.000 A*S (LA40 HP 8000 N)</p> <p>LA44: 10.000 cycles in life test equals: 5.000.000 A*S</p> <p>BL1: 10.000 cycles in life test equals: 2.200.000 A*S</p> <p>LC1: 10.000 cycles in life test equals: 3.290.000 A*S</p> <p>LC3: 10.000 cycles in life test equals: 6.400.000 A*S</p> <p>The work indicator on each actuator can be reset by pressing </p>
<p>Overloads</p> <p></p> <p>0</p>	<p>Total number of overloads reached inwards and outwards. </p> <p>If more actuators are moving at the same time during an overload situation the overload figure will be counting on all running actuators</p>
<p>Acc.Duration</p> <p>Σ[s]</p> <p>85</p>	<p>Accumulated activity duration time in seconds.</p>



Symbol	Explanation
<p>IN</p>  <p>8</p>	Read out how many times a control button has been activated inwards
<p>OUT</p>  <p>7</p>	Read out how many times a control button has been activated outwards
	Numbers of errors per actuator. The indicator can only be reset by exchanging the actuator through SDT2 "Actuator" menu.
 <p>No. of resets</p> <p>0</p>	<p>Activate this button for update of data when the actuator has been replaced. When a new actuator has been connected it can be re-initiated to start from zero. If "replace" is not activated the COXX believes the first actuator is still connected and the wrong data will be shown.</p> <p>Please notice: A new actuator connected is only visible in the Service Data Tool view if it has been running in outward or inward direction – A*S>0.</p> <p>When you activate this button it will also reset the number of errors per actuator – see the summary section.</p> <p>The number "0" indicates the number of resets per actuator. Before resetting is complete there will be dialogue boxes like...</p> <div data-bbox="395 974 756 1547" style="border: 1px solid black; padding: 10px;">  </div>

How to conclude on the service information


Total work:

Please contact the manufacturer of the bed application or medical application e.g. couch / table or chair for treatment and examination. The manufacturer will decide when it is appropriate to consider exchanging the actuator.

Overload:

If overload has occurred it is recommended to consider stronger lifting equipment with higher working load for the particular patient / institution.


Controls view



Symbol	Explanation
A Openbus "H-Frame" Input 	OpenBus Signals
B Openbus "A-Frame" Led's 	OpenBus Diode Signals
C Openbus "A-Frame" Input 	OpenBus Button Signals
D ID4-ID7 Service (ID4) ExBits (ID5) Res (ID6)	Additional OpenBus Signals overview



Report view




Example given with a 4-actuator system

Symbol	Explanation
<p>A</p> <p>System P.O <input type="text" value="12345678"/></p> <p>System Item No. <input type="text" value="87654321"/></p> <p>System S.W <input type="text" value="13579"/></p> <p>Ch1 Item No. <input type="text" value="LA23-1"/></p> <p>Ch2 Item No. <input type="text" value="LA23-2"/></p> <p>Ch3 Item No. <input type="text" value="LA23-3"/></p> <p>Ch4 Item No. <input type="text" value="LA23-4"/></p> <p>Item No. <input type="text" value="MJB500-1023"/></p> <p>Item No. <input type="text"/></p> <p>Item No. <input type="text"/></p> <p>Item No. <input type="text"/></p> <p>Item No. <input type="text"/></p> <p>Software Number <input type="text" value="SW20100290V0-126602856"/></p>	Write the information which is on the label of the product.
<p>B</p> <p>Inspected By <input type="text" value="Service Technician"/></p> <p>Company <input type="text" value="Service Company XX"/></p>	Fill in the information; Inspected by, Company
<p>C</p> <p>Action Taken <input type="text" value="This text will only be saved in the document – pdf or csv file."/></p>	Description made in this box will be read out when the report is saved.
<p>D</p> <p>Notes: <input type="text" value="003: Actuator 2 replaced"/> <input type="button" value="v"/></p> <p>This Service Visit <input type="text" value="007: Battery replaced"/> <input type="button" value="v"/></p> <p><input type="text"/> <input type="button" value="v"/></p> <p><input type="text"/> <input type="button" value="v"/></p> <p style="text-align: center;">x</p> <p style="font-size: small;">Service code has been saved successfully to the control box.</p> <p style="text-align: center;"><input type="button" value="OK"/></p>	Choose up to 4 notes after the service visit. Messages for the next service visit - will be stored and readable at next visit. Saving is confirmed.
<p>E</p> <p>Notes: Previous Service Visit <input type="text" value="This information will be saved in the control box until the next report is saved."/></p>	Up to 4 notes readable from the previous service visit.




Example given with a 4-actuator system



F



G




When the application has been checked and the information completed, press “save” to save as PDF, CSV, send to Email. Relevant actuator statistics and system data will be included in the report.

Please notice that outlook has to be available. Outlook will open automatically when report download is started.





If the COXX is connected to a Battery, Lithium Ion battery – all battery data will also be saved



Advanced settings view




Example given with a 6-actuator system

Symbol	Explanation
<p>A</p> <p>ACTUATOR STATISTICS</p> 	<p>Remember to update actuator info if the COXX is replaced. In this way you maintain the statistical service data information on the actuators.</p>
<p>B</p> <p>SYSTEM STATISTICS</p> 	<p>Read out data related to unique customised functions</p>
<p>C</p> 	<p>Get/Retrieve data from CB: When entering new data the data will be visible with bold script. "Get Data from CB" undo changes and retrieve existing data from the control box.</p>
<p>D</p> <p>Please remember to save report with existing data before exchanging the CO61</p> 	<p>Remember to save report with existing data before exchanging the Control box. Save new settings: This will reset the above information with the new data filled in.</p>



Troubleshooting view




Example given with a 6-actuator system

Symbol	Explanation																																																																																											
<p>A</p> <ul style="list-style-type: none"> <input type="radio"/> Control box error 0 <input type="radio"/> Controls error 2 <input type="radio"/> Duty cycle guard 0 	<p>The LED is orange if there is a failure or error. The orange triangle symbol is shown next to the channel where the failure is detected.</p> <p>The total number of errors is stated and can only reset by replacing the control box and press the following button in the Advanced settings section</p> <div style="border: 1px solid black; padding: 5px; width: fit-content; margin: 10px auto;">Get Data from CB</div>																																																																																											
<p>B</p> <ul style="list-style-type: none"> <input type="radio"/> Actuator error 	<p>The LED is orange if there is an actuator error. Please see the actuator statistics for the total number of errors per actuator.</p>																																																																																											
<p>C</p> <table border="1" style="width: 100%; border-collapse: collapse; font-size: 8px;"> <thead> <tr> <th colspan="7">ACTUATOR STATISTICS</th> </tr> <tr> <th></th> <th>1</th> <th>2</th> <th>3</th> <th>4</th> <th>5</th> <th>6</th> </tr> </thead> <tbody> <tr> <td>Cutoff Limit OUT</td> <td>3</td> <td>3</td> <td>3</td> <td>3</td> <td>9</td> <td>9</td> </tr> <tr> <td>Cutoff Limit IN</td> <td>3</td> <td>3</td> <td>3</td> <td>3</td> <td>9</td> <td>9</td> </tr> <tr> <td>Total Overloads</td> <td>0</td> <td>0</td> <td>0</td> <td>0</td> <td>0</td> <td>0</td> </tr> <tr> <td>Total Work</td> <td>51</td> <td>29</td> <td>17</td> <td>22</td> <td>1235</td> <td>1204</td> </tr> <tr> <td>Total Acc. Duration</td> <td>229</td> <td>127</td> <td>77</td> <td>96</td> <td>1271</td> <td>1292</td> </tr> <tr> <td>Activations OUT</td> <td>73</td> <td>32</td> <td>11</td> <td>32</td> <td>162</td> <td>197</td> </tr> <tr> <td>Activations IN</td> <td>57</td> <td>34</td> <td>22</td> <td>20</td> <td>143</td> <td>151</td> </tr> <tr> <td>Actuator error type 1</td> <td>0</td> <td>0</td> <td>0</td> <td>0</td> <td>3</td> <td>5</td> </tr> <tr> <td>Actuator error type 2</td> <td>0</td> <td>0</td> <td>0</td> <td>0</td> <td>1</td> <td>1</td> </tr> <tr> <td>Actuator error type 3</td> <td>0</td> <td>0</td> <td>0</td> <td>0</td> <td>0</td> <td>0</td> </tr> <tr> <td>Position Lost</td> <td>0</td> <td>0</td> <td>0</td> <td>0</td> <td>0</td> <td>0</td> </tr> </tbody> </table>	ACTUATOR STATISTICS								1	2	3	4	5	6	Cutoff Limit OUT	3	3	3	3	9	9	Cutoff Limit IN	3	3	3	3	9	9	Total Overloads	0	0	0	0	0	0	Total Work	51	29	17	22	1235	1204	Total Acc. Duration	229	127	77	96	1271	1292	Activations OUT	73	32	11	32	162	197	Activations IN	57	34	22	20	143	151	Actuator error type 1	0	0	0	0	3	5	Actuator error type 2	0	0	0	0	1	1	Actuator error type 3	0	0	0	0	0	0	Position Lost	0	0	0	0	0	0	<p>Overview of the number of errors per channel. </p> <p>When an actuator has been replaced and the data are updated by pressing the refresh button in the actuator section – the number of failure or error per channel will also be updated</p>
ACTUATOR STATISTICS																																																																																												
	1	2	3	4	5	6																																																																																						
Cutoff Limit OUT	3	3	3	3	9	9																																																																																						
Cutoff Limit IN	3	3	3	3	9	9																																																																																						
Total Overloads	0	0	0	0	0	0																																																																																						
Total Work	51	29	17	22	1235	1204																																																																																						
Total Acc. Duration	229	127	77	96	1271	1292																																																																																						
Activations OUT	73	32	11	32	162	197																																																																																						
Activations IN	57	34	22	20	143	151																																																																																						
Actuator error type 1	0	0	0	0	3	5																																																																																						
Actuator error type 2	0	0	0	0	1	1																																																																																						
Actuator error type 3	0	0	0	0	0	0																																																																																						
Position Lost	0	0	0	0	0	0																																																																																						
<p>D</p> <div style="border: 1px solid black; padding: 5px; margin-bottom: 10px; text-align: center;">Download and send log file</div> <div style="border: 1px solid black; padding: 10px; text-align: center; margin-bottom: 10px;">Download is in progress. Please wait...</div> <div style="border: 1px solid black; padding: 10px; text-align: center;">The report has been sent to email.</div> <div style="text-align: right; margin-top: 5px;"><input type="button" value="Ok"/></div>	<p>Press this button to download and send a log file with all actuator statistics and troubleshooting data.</p> <p>The log file will be sent to your outlook for easy forward to your LINAK contact for further support</p> <p>Please notice that outlook has to be available. Outlook will open automatically when report download is started.</p> <p>Status information and confirmation of successful download of Troubleshooting report will be shown.</p>																																																																																											




PJ2 Junction box view






	Symbol	Explanation
	<i>Please note that this section is only visible when PJ2 is connected to a LINAK control box. Data are slowly updated.</i>	
A	PJ2 details	Automatic read-out for the PJ2 for easy identification.






Battery view



Symbol	Explanation
	<i>Please note that this section is only visible when COXX is connected to a LINAK battery.</i>
A	 <p>This symbol indicates that COXX is connected to mains power.</p>
	 <p>This symbol indicates that COXX is operating on a Lithium Ion battery. Battery Status is indicated by % and colour: Green: Battery fully charged (approx. 100-40% capacity remaining) Yellow: Charging recommended (approx. 40-21% capacity remaining) Orange: Low critical battery level. (Depending on the COXX SW there will be audio signal when the hand control is activated. Limited the actuator function).</p>
B	<p>Battery details</p> <p>Automatic read-out of the battery details for easy identification of the battery</p>
C	<p>Battery life</p>  <p>Battery life is indicated by % and colour: Green: Excellent battery life (approx. 100-70% capacity remaining) Yellow: Good battery life (approx. 69-30% capacity remaining) Orange: Critical battery life - replace battery (below 30% capacity remaining)</p>
D	<p>Battery health</p> <p>The use of the battery affects the battery health and the number of incidents of the following can also help explaining the battery health and the length of the battery life. These statistical data are also useful for service evaluation.</p>
E	<p>Battery use</p> <p>This section indicates the use of the battery.</p> <p>The work indicator for the battery measures via ampere usage*seconds in use. The work indicator gives a very good indication of how much the battery is worn.</p> <p>These statistical data are also useful for service evaluation. For BA21 battery, 500 full charging/decharging cycles is approx. 2.700.000 A *S. However, this is only a rough guideline as a lot of factors will affect the BA21 lifetime, e.g. how the battery is used and charged.</p> <p>See the Battery Section for further details.</p>


BA22 Battery Service view





Symbol	Explanation
	<p>Please notice! This Service Data Tool section is only to be used when support of LINAK BA22 lithium Ion battery. Please notice that outlook has to be available. Outlook will open automatically when report download is started. If you have a CO-Link system, it is possible to read-out battery statistics from two batteries connected to the primary COXX and the secondary COXX respectively.</p> <p>Choose the BA22 view Battery for read-out of statistics connected to the COXX primary control box.</p> <p>Choose the BA22 view Battery secondary for read-out of statistics connected to the COXX secondary control box.</p> <p>Having a system where the battery is connected to a PJ2, the duration for generating the report can take up to 10 min.</p>
A	<p>Press "save" to download a battery log file as CSV and send to Email. The Battery details and relevant Battery statistics will be included in the report. These data are useful for further technical analysis of the battery.</p>  <p>Download status indication.</p> 
B	<p>Battery details</p> <p>Automatically read-out of the battery details for easy identification of the battery</p>
















BA16/19 Battery Service view



Symbol	Explanation
	<p>Please notice! This Service Data Tool section is only to be used when support of LINAK BA16 / BA19 Lead acid batteries.</p> <p>Please notice that outlook has to be available. Outlook will open automatically when report download is started. If you have a CO-Link system, the Service Data Tool it is possible to read-out battery statistics from two batteries connected to the primary COXX and the secondary COXX respectively.</p> <p>Choose the BA16/BA19 view for read-out of statistics from BA16/BA19 connected to the COXX primary control box.</p> <p>Choose the BA16/BA19 2nd view for read-out of statistics from BA16/BA19 connected to the COXX secondary control box.</p>
A	<p>Press "save" to download a battery log file as CSV and send to Email. The Battery details and relevant Battery statistics will be included in the report. These data are useful for further technical analysis of the battery.</p>  <p>Download status indication.</p> 
B	<p>Battery details</p> <p>Automatically read-out of the battery details for easy identification of the battery</p>



How to change the application picture, company logo and information PDF-file


	Symbol	Explanation
		Click on the application picture or company logo with the right mouse button to select another picture or logo. Please use the file type JPEG or PNG.
		Choose the application picture from your desktop. The size of the new application picture or company logo should be: Application picture: Width 138 pixels, Height 112 pixels Company logo: Width 150 pixels, Height 18 pixels
1		To change the pdf files: help or info –  Open Computer 
2		Open the drive including "LINAK"
3		Open SDT2
4		Open Devices
5		Open the software-map for COXX
6		Click on devices, COXX and views
7		Click on PDF to change a pdf file in one of the sections or to add an info file
8		Rename the existing pdf file. For example from "Info site" to "OLD" Remember to save or write the name you just have changed. In this example you have to remember "Info site"
9		Copy the file you want to add to the folder and rename the new pdf with the name of the pdf you want to replace. If we take the example from point 8, you have to change it to "Info site"

SDT2 for CO-Link systems

System overview for CO-Link systems

The LINAK IB3 tool (the SDT “box”) also allows for service data read-out from LINAK CO-Link systems with a modular (also called RJ50) connection. Below is a diagram of how this IB3 tool is connected to a LINAK CO-Link system via a Modular Junction Box (MJB). For CO-Link, it is not possible to connect the yellow connector directly to the CO control box.


An MJB is required to maintain the CO-Link communication, and an MJB is recommended so other components can still be connected and ready to send commands, both to wake up the system for the SDT2 program, but also to run the system for troubleshooting purposes.




CO-Link system explanation

When Service Data Tool is opened there will be access to the versions below. COXX, BA16/BA19 and BA22 are the primary views. Secondary views are only used for CO-Link systems.

Example:



It is easy to identify the selected view by looking in the upper left corner – here it is the COXX primary view used for the CO-Link system.



	Symbol	Explanation
A		<p>CO-Link is used when two control boxes e.g. COXX are connected in a system with control of up to 8 movements / actuators.</p> <p>Look into the device information to see if the service data are read out from the COXX primary or secondary.</p> <p>The COXX primary and COXX secondary views have exact the same layout and information available. Only difference is that data are read-out from two different control boxes – COXX.</p>
B		<p>Service Data Tool COXX primary shows service data for up to 4 actuators.</p> <p>Service Data Tool for COXX secondary shows service data for up to 4 actuators.</p> <p>Service Data Tool can only be opened for the COXX primary OR the COXX secondary – one at the time.</p> <p>Choose between the two Service Data Tool versions</p> <p>COXX (change) – Service Data Tool for COXX Primary</p> <p>COXX Secondary (change) – Service Data Tool for COXX Secondary</p>
C	Status bar	<p>Service Data Tool is either IDLE or reading data.</p> <p>The message means that SDT is reading data – please wait</p> <p>Reading Extended openbus device 8 register20</p>




SDT2 for LIFT50 systems

System overview for LIFT50

The LINAK IB3 tool (the SDT “box”) can also be used for service data read-out from the LIFT50 control box with a modular (also called RJ50) connection. Below is a diagram of how this IB3 tool is connected to a standard LIFT50 system via a Modular Junction Box (MJB).

While it is also possible to connect the yellow connector directly to the CO control box, it is advisable to do so through an MJB so other components can still be connected and ready to send commands, both to wake up the system for the SDT2 program, but also to run the system for troubleshooting purposes.



Information


Please note that the following view is now called “LIFT” instead of “COL50” in SDT ver. 3 and above and is an included view (does not need to be imported first).









Reading out service data on a PC


Ensure that the COL50 or battery BAL50 view is initiated by pressing the menu shown.

Please contact your local LINAK Supplier for support if this is not the case.




LIFT50 Service Data Tool sections









	Symbol	Explanation
A		For quick and easy overview of actuators info, service counters and battery state. Please notice: A new actuator connected is only visible in the Service Data Tool view if it has been running in outward or inward direction – $A*S > 0$.
B		For detailed information about the actuators' statistical service data. For refreshing data from the control unit connected. For update of information when replacing the actuator.
C		For detailed information about hand control signals and codes.
D		For service reporting, production number, software number, item number of the application. Saving complete and relevant information for actuators data and service.
E		Intended for trained and authorised service technicians only. For update of actuator info if the control is replaced to maintain service data. For update of service settings and cutoff limit settings.
F		For information about battery life, charging state and use.
G	Status bar	Status is either IDLE or Reading Extended OpenBus device 8 register 20 It means that SDT is reading data – please wait.



Sections are updated every second to keep data updated. Please re-start SDT if the control box has powered down or disconnected. For further help and recommendations on each section, please press the help icon 

Please re-start Service Data Tool if hardware has been exchanged or dismantled during monitoring.


Summary view



Symbol	Explanation
A	 <p>This symbol indicates that COL50 is connected to mains power and the current battery status.</p>
	 <p>This symbol indicates that COL50 is operating on a BAL50 Lithium Ion battery.</p>
B	<p>Total number of service counter reset: 0</p> <p>Amount of service counters re-set in total during application lifetime.</p>
C	 <p>Graphical overview of actual service counter setting and state. Service counter is shown in either days or cycles. One cycle is defined as: 5 sec. outwards, 2 sec. inwards operation and current consumption > 2 A. Service counter state is indicated by colours. when service counter settings are getting expired, indicator moves into yellow and red area.</p>
D	<p>Selected or de-selected the actuators you want to monitor</p>
	 <p>The actuator position is shown with numbers.</p>
	 <p>The green arrow indicates that the actuator is working and shows which direction, which is also written. Status can be; Direction In, Direction Out, EOS in, EOS out or IDLE</p>
	 <p>If the actuator has lost its position, POS LOST indicator will occur in the status bar</p>
E	 <p>Help manual of what the different functions means.</p>
F	 <p>The customer can add their own information as PDF-file.</p>

G		A picture of the control box
H		The customer can change the application picture.
I	Company logo	The customer can add their own LOGO.




Actuator view



Example given with a 3-actuator system

	Symbol	Explanation
A		<p>To ensure that the application will stop if the current draw exceeds the preset limits.</p> <p>These values are factory default setting.</p>
		<p>Total work on the actuator (A*S): Work indicator for the actuators measures via ampere usage *seconds in use. The work indicator gives a very good indication of how much the actuator is worn.</p> <p>LA20: 10.000 cycles in life test equals: 1.300.000 A*S LA23: 10.000 cycles in life test equals: 1.400.000 A*S LA27: 10.000 cycles in life test equals: 3.700.000 A*S LA31: 10.000 cycles in life test equals: 3.000.000 A*S LA34: 10.000 cycles in life test equals: 4.300.000 A*S LA43: 10.000 cycles in life test equals: 2.500.000 A*S LA40: 10.000 cycles in life test equals: 3.800.000 A*S (LA40 6000N) LA40: 10.000 cycles in life test equals: 4.400.000 A*S (LA40 PL 8000N) LA40: 10.000 cycles in life test equals: 2.900.000 A*S (LA40 HP 8000N) LA44: 10.000 cycles in life test equals: 5.000.000 A*S BL1 : 10.000 cycles in life test equals: 2.200.000 A*S LC3 : 10.000 cycles in life test equals: 6.400.000 A*S</p> <p>The work indicator on each actuator can be reset by pressing </p>
		<p>Total number of overloads reached inwards and outwards. </p> <p>If more actuators are moving at the same time during an overload situation the overload figure will be counting on all running actuators</p>



Symbol	Explanation
 $\Sigma[s]$ 85	Accumulated activity duration time in seconds.
IN  8	Read out how many times a control button has been activated inwards
OUT  7	Read out how many times a control button has been activated outwards

How to conclude on the service information


Total work:




Please contact the manufacturer of the bed application or medical application e.g. couch / table or chair for treatment and examination. The manufacturer will decide when it is appropriate to consider exchanging the actuator.

Overload:

If overload has occurred it is recommended to consider stronger lifting equipment with higher working load for the particular patient / institution.


COL50 Controls view



	Symbol	Explanation
A	Openbus "H-Frame" Input 	OpenBus Signals
B	Openbus "A-Frame" Led's 	OpenBus Diode Signals
C	Openbus "A-Frame" Input 	OpenBus Button Signals
D	ID4-ID7 Service (ID4) ExBits (ID5) Res (ID6)	Additional OpenBus Signals overview



COL50 Report view



Example given with a 3-actuator system


Symbol	Explanation
<p>A</p> <p>System P.O. <input type="text"/></p> <p>System Item No. <input type="text"/></p> <p>System S.W. <input type="text"/></p> <p>Controlbox serial No. <input type="text" value="1234567890"/></p> <p>Ch1 Item No. <input type="text" value="LA40"/></p> <p>Ch2 Item No. <input type="text" value="LA23"/></p> <p>Ch3 Item No. <input type="text" value="LA20"/></p> <p>Item No. <input type="text" value="J12345"/></p> <p>Item No. <input type="text" value="J23456"/></p> <p>Item No. <input type="text"/></p> <p>Item No. <input type="text"/></p> <p>Item No. <input type="text"/></p> <p>Software Number <input type="text" value="SW19700120V0-212844722"/></p>	<p>Write the information which is on the label of the product.</p>
<p>B</p> <p>Inspected By <input type="text" value="Service Technician"/></p> <p>Company <input type="text" value="Service Company XX"/></p>	<p>Fill in the information; Inspected by, Company</p>
<p>C</p> <p>Action Taken <input type="text" value="This text will only be saved in the document – pdf or csv file."/></p>	<p>Description made in this box will be read out when the report is saved.</p>
<p>D</p> <p>Notes:</p> <p>This Service Visit</p> <p><input type="text" value="003: Actuator 2 replaced"/> <input checked="" type="checkbox"/></p> <p><input type="text" value="007: Battery replaced"/> <input checked="" type="checkbox"/></p> <p><input type="text"/> <input checked="" type="checkbox"/></p> <p><input type="text"/> <input checked="" type="checkbox"/></p>	<p>Choose up to 4 notes after the service visit.</p> <p>Messages for the next service visit - will be stored and readable at next visit. Saving is confirmed.</p> <p style="text-align: center;">×</p> <p style="text-align: center;">Service code has been saved successfully to the control box.</p> <p style="text-align: center;"><input type="button" value="OK"/></p>



	Symbol	Explanation
E	<p>Notes: Previous Service Visit</p> <div data-bbox="354 174 584 271" style="border: 1px solid black; padding: 2px; width: fit-content;"> This information will be saved in the control box until the next report is saved. </div>	<p>Up to 4 notes readable from the previous service visit.</p>
F	<p>Service periode</p> <div data-bbox="411 327 691 358" style="border: 1px solid black; padding: 2px; width: fit-content;"> 360 days </div> <div data-bbox="411 371 691 403" style="border: 1px solid black; padding: 2px; width: fit-content;"> 8000 cycles </div>	<p>Defined service periode in days and cycles. One cycle is defined as: 5 sec. outwards and 2 sec. inwards operation. Service counter state is indicated by colours. When service counter settings are getting expired, indicator moves into yellow and red area.</p>




Example given with a 3-actuator system



	Symbol	Explanation
G		When the application has been checked and the information completed, press “save” to save as PDF, CSV, send to Email. Relevant actuator statistics and system data will be included in the report.
H		<p>Please notice that outlook has to be available. Outlook will open automatically when report download is started.</p> <p>If the COL50 is connected to BAL50 battery - all battery data will be saved.</p>






COL50 Advanced settings



Symbol	Explanation																																																
<p>A ACTUATOR STATISTICS</p> <table border="1"> <tr> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>Cutoff Limit OUT</td> <td>↑ 8</td> <td>5,5</td> <td>5,5</td> </tr> <tr> <td>Cutoff Limit IN</td> <td>↓ 8</td> <td>3</td> <td>3</td> </tr> <tr> <td>Overloads OUT</td> <td>↑ 0</td> <td>0</td> <td>0</td> </tr> <tr> <td>Overloads IN</td> <td>↓ 0</td> <td>0</td> <td>0</td> </tr> <tr> <td>Work OUT</td> <td>↑ 91</td> <td>13</td> <td>3</td> </tr> <tr> <td>Work IN</td> <td>↓ 85</td> <td>33</td> <td>13</td> </tr> <tr> <td>Acc. Duration OUT</td> <td>↑ 55</td> <td>12</td> <td>4</td> </tr> <tr> <td>Acc. Duration IN</td> <td>↓ 44</td> <td>22</td> <td>8</td> </tr> <tr> <td>Activations OUT</td> <td>↑ 6</td> <td>4</td> <td>1</td> </tr> <tr> <td>Activations IN</td> <td>↓ 10</td> <td>3</td> <td>2</td> </tr> <tr> <td>Total Cycles</td> <td>↑ 4</td> <td></td> <td></td> </tr> </table>					Cutoff Limit OUT	↑ 8	5,5	5,5	Cutoff Limit IN	↓ 8	3	3	Overloads OUT	↑ 0	0	0	Overloads IN	↓ 0	0	0	Work OUT	↑ 91	13	3	Work IN	↓ 85	33	13	Acc. Duration OUT	↑ 55	12	4	Acc. Duration IN	↓ 44	22	8	Activations OUT	↑ 6	4	1	Activations IN	↓ 10	3	2	Total Cycles	↑ 4			<p>Remember to update actuator info if the COL50 is replaced. In this way you maintain the statistical service data information on the actuators.</p> <p>Cutoff Out and Cutoff In values are default factory settings for current Cutoff limits.</p> <p>Valid for CH1 in case of O-values in learn mode fields.</p>
Cutoff Limit OUT	↑ 8	5,5	5,5																																														
Cutoff Limit IN	↓ 8	3	3																																														
Overloads OUT	↑ 0	0	0																																														
Overloads IN	↓ 0	0	0																																														
Work OUT	↑ 91	13	3																																														
Work IN	↓ 85	33	13																																														
Acc. Duration OUT	↑ 55	12	4																																														
Acc. Duration IN	↓ 44	22	8																																														
Activations OUT	↑ 6	4	1																																														
Activations IN	↓ 10	3	2																																														
Total Cycles	↑ 4																																																
<p>B</p> <table border="1"> <tr> <td>Service settings</td> <td></td> </tr> <tr> <td>Days</td> <td>360</td> </tr> <tr> <td>Cycles</td> <td>8000</td> </tr> </table>	Service settings		Days	360	Cycles	8000	<p>Service interval current settings for this Lift system. When entering new data, the data will be visible in bold script. Press save data to control box to update.</p>																																										
Service settings																																																	
Days	360																																																
Cycles	8000																																																
<p>C</p> <table border="1"> <tr> <td>Simple Learn</td> <td></td> </tr> <tr> <td>Cutoff Limit OUT</td> <td>1,7</td> </tr> <tr> <td>Advanced Learn</td> <td></td> </tr> <tr> <td>Cutoff Limit 1</td> <td>0</td> </tr> <tr> <td>Cutoff Limit 2</td> <td>0</td> </tr> <tr> <td>Cutoff Limit 3</td> <td>0</td> </tr> <tr> <td>Cutoff Limit 4</td> <td>0</td> </tr> <tr> <td>Cutoff Limit 5</td> <td>0</td> </tr> <tr> <td>Cutoff Limit 6</td> <td>0</td> </tr> <tr> <td>Cutoff Limit 7</td> <td>0</td> </tr> <tr> <td>Cutoff Limit 8</td> <td>0</td> </tr> </table>	Simple Learn		Cutoff Limit OUT	1,7	Advanced Learn		Cutoff Limit 1	0	Cutoff Limit 2	0	Cutoff Limit 3	0	Cutoff Limit 4	0	Cutoff Limit 5	0	Cutoff Limit 6	0	Cutoff Limit 7	0	Cutoff Limit 8	0	<p>CH1 current cutoff limit settings for this lift system. Data can either be a result of manual entries or be created via learn mode function on your handset.</p> <p>Press save after learn.</p> <p>If learn mode is mandatory for operation, values are to be filled.</p> <p>Actuators with feedback required for Advanced Lean.</p>																										
Simple Learn																																																	
Cutoff Limit OUT	1,7																																																
Advanced Learn																																																	
Cutoff Limit 1	0																																																
Cutoff Limit 2	0																																																
Cutoff Limit 3	0																																																
Cutoff Limit 4	0																																																
Cutoff Limit 5	0																																																
Cutoff Limit 6	0																																																
Cutoff Limit 7	0																																																
Cutoff Limit 8	0																																																



	Symbol	Explanation
D	 FOR EXPERTS ONLY	<p>Export function to upload learn mode values into software configurator.</p> <p>The file will be save to PC desktop.</p> <p>File to be send to LINAK representativ for general settings.</p>
E		<p>Get/Retrieve data from CB: When entering new data the data will be visible with bold script. "Get Data from CB" undo changes and retrieve existing data from the control box.</p>
F	<p>Please remember to save report with existing data before exchanging the CO61</p> 	<p>Remember to save report with existing data before exchanging the Control box. Save new settings: This will reset the above information with the new data filled in.</p>



COL50 Battery View

LINAK Battery v1.7



ACTUAL CHARGING STATE
 1.5 Ah 73 %
 On Mains

BATTERY HEALTH
 Battery health: Excellent
 Number of under voltage: 0
 Number of discharge overcurrent: 0
 Number of undertemperature: 0
 Number of overtemperature: 0

BATTERY DETAILS
 Battery capacity designed: 2.1 Ah
 SW number / version: SW00999503V0-52865
 Production week / year: 25/45/535
 PO number and series number: 4294967295-65535


BATTERY USE
 Actual temp. in battery box: 27°C 80°F
 Actual pack current: 12.2 A
 Total A/S charged into battery: 15144 A/S



Company Logo

Symbol	Explanation
<p>A</p>  	<p>This symbol indicates that COL50 is connected to mains power and the current battery status.</p> <p>Battery Status is indicated by % and Ah. If battery capacity is less than 2 cycles the COL50 will hoot.</p>
B Battery details	Automatically read-out of the battery details for easy identification of the battery
C Battery health	The use of the battery affects the battery health and the number of incidents of the following can also help explaining the battery health and the length of the battery life. These statistical data are also useful for service evaluation.
D Battery Use	<p>This section indicates the use of the battery.</p> <p>The work indicator for the battery measures via ampere usage *seconds in use. The work indicator gives a very good indication of how much the battery is worn.</p> <p>These statistical data are also useful for service evaluation. For BAL50 battery 500 full charging /decharging cycles is approx. 3.510.000 A x S. However, this is only a rough guideline as a lot of factors will affect the Lift battery life time e.g. how the battery is used and charged. See the Battery Section for further details.</p>
















BAL50 Battery Service View



Symbol	Explanation
	Please notice that outlook has to be available. Outlook will open automatically when report download is started.
A	<p>Press "save" to download a battery log file as CSV and send to Email. The Battery details and relevant Battery statistics will be included in the report. These data are useful for further technical analysis of the battery.</p>  <p>Download status indication.</p> 
B	Battery details
	Automatically read-out of the battery details for easy identification of the battery



How to change the application picture, company logo and information PDF-file

	Symbol	Explanation
		Click on the application picture or company logo with the right mouse button to select another picture or logo. Please use the file type JPEG or PNG.
		Choose the application picture from your desktop. The size of the new application picture or company logo should be: Application picture: Width 138 pixels, Height 112 pixels Company logo: Width 150 pixels, Height 18 pixels
1		To change the pdf files: help or info –   Open Computer
2		Open the drive including "LINAK"
3		Open SDT2
4		Open Devices
5		Open the software-map for COL50
6		Click on views, COL50
7		Click on PDF to change a pdf file in one of the sections or to add an info file
8		Rename the existing pdf file. For example from "Info site" to "OLD" Remember to save or write the name you just have changed. In this example you have to remember "Info site"
9		Copy the file you want to add to the folder and rename the new pdf with the name of the pdf you want to replace. If we take the example from point 8, you have to change it to "Info site"


SDT2 for JUMBO Care systems

System overview for JUMBO Care

The LINAK IB3 tool (the SDT "box") can also be used for service data read-out from the JUMBO Care control box with a DIN connection. Below is a diagram of how this IB3 tool is connected to a standard LIFT50 system via a Modular Junction Box (MJB).

While it is also possible to connect the DIN connector directly to the CO control box, it is advisable to do so using the T-cable (0015412-A) so that other components can still be connected and ready to send commands, both to wake up the system for the SDT2 program, but also to run the system for troubleshooting purposes.

Service Data Tool and JUMBO Care 2nd generation, 3 channel solution:



Information

Please note that the Service Data Tool version 2.5.0 or newer is needed. The JUMBO Care article number should begin with CBJxxxxxxxxxxxx.

Please remember to order cables and the IB300001 box for data read-out from JUMBO Care to laptop.

The T-Cable 0015412-A is needed if you do not have an activation button on the JUMBO Care control box.



Article number	Product
IB300001	Box for data read out from JUMBO Care to laptop
0964987-A	Cable for usage of Service Data Tool (from JUMBO Care to IB300001 box)
0015412-A	T-Cable for usage of hand control and Service Data Tool. (from JUMBO Care to hand control and IB300001 box)
0964549-C	T-Cable, Stereo Jack for CH1 to enable 3 channel solution



Service intervals

What is the service interval presetting?

A standard JUMBO Care control box will be pre-set with a service interval of: 12 months / 8000 cycles, whichever comes first.

How do I see on the control box that it is time for service?

The notice about service need depends upon the JUMBO Care version:

- No indicators on the front cover: one single beep telling that it is time for service
- Diodes on the front cover: The service diode will light up and one single beep will tell that it is time for service
- Display: The display will show the service symbol and one single beep will tell that it is time for service.

When does the service time start counting?

The service time will start to count down from the day the control box is produced.


How to change to another service interval?

The service interval can be changed via a special LINAK handset (Item number HB8646V2010 + 71). To change the service interval you simply plug in the handset in the JUMBO Care and press the button. You can choose from a number of months or cycles between services. You can also choose 'NO months' or 'NO cycles'.

E.g. Push M = 24, Push \updownarrow 10000.

The JUMBO Care is now set to indicate service need every 24 month or at 10000 cycles, whichever comes first.

After having set a different service interval, the control box will provide an audio sound (as a receipt).



Read out service data on a JUMBO Care with display

When you have a JUMBO Care with display, it is possible to have some basic service data on the display. To get this information on the display, press the “lifting arm up” button on your LINAK hand control or control box (short press ½ second).


The information that appears on the display is;







CBJ Care read out service data



HB80




HB70

	12034	→ Total cycles made by the lifting actuator (channel 1)
	1257000	→ Total work made by the lifting actuator (channel 1)
	7	→ Total number of overloads (channel 1)
	90/360	→ Days since last service/days between services

Reading out service data on a PC

Ensure that the JUMBO Care view is initiated by pressing the menu shown.

Please contact your local LINAK supplier for support if this is not the case.



Information

Please note that the "Jumbo Care_2gen" view is now simply called "JumboCare" for both 1st and 2nd generation as of SDT2 ver. 3 and above. If you need this view, please contact your local LINAK subsidiary.

The Service Data Tool for Jumbo Care is divided into sections:

Section	Used for.....
Summary	For quick and easy overview of service indicator, battery status and actuator statistics
Actuators	For detailed information about the actuators' statistical service data
Control box	For detailed information about the control box and battery state
Controls_InputOutput	For detailed information about hand control signals and codes
Report	For service reporting and resetting of service counters via hand control after saving report
Advanced settings	Intended for trained and authorised service technicians only. For change of service settings and update of actuator info if the control box is replaced to maintain service data. For update of information when replacing the actuator

For further help and recommendations on each section, please press the help icon or see next page.







Actuator view

The screenshot shows the 'LINAK JUMBO CARE Actuators view' with the following data:

Actuator	Status	Cycles	Work (A*S)	Overloads In	Overloads Out	Acc. Duration In	Acc. Duration Out
1	Idle	2	270	2	3	29	32
2	Idle	4	0	0	0	1	0
3	Idle	5	0	0	0	0	3

Additional data from the screenshot:

- Connection Info: CDM4 (v100.1.7.1), OpenbusDK
- Devices: CBUC, Device: 4, SPT: 7, app. SW: 1000003, Version: 0.9800
- Actuator 1 Cutoff limit: In 10,7, Out 10,7
- Actuator 2 Cutoff limit: In 4,0, Out 4,0
- Actuator 3 Cutoff limit: In 4,0, Out 4,0
- Reading Extended openbus device 4 register76
- Last Updated 13:49:20

Symbol	Service data description and recommendations
 1  Status: Idle	<p>If the overload symbol is shown next to the actuator there is currently an overload situation.</p> <p>Reduce the load on the actuator. If a buffer actuator is used, the amount of overloads is an expression of how many times the actuator has run into end-stop and been overloaded.</p>
Overload  4 Shown in the summary section	<p>Total overloads CH1:</p> <p>Counts the number of times the actuator on channel 1 has been overloaded. If the actuator has been overloaded, we recommend lifting equipment with higher load capacity. The indicator can only be reset by exchanging the actuator through Service Data Tool "actuator" menu.</p> <p>Channels 2 and 3: If a buffer actuator is used, the amount of overloads is an expression of how many times the actuator has run into end-stop and been overloaded</p>
Work A * S 481	<p>Total work CH1 (A*S):</p> <p>Work indicator for the actuator on channel 1 measures via ampere usage * seconds in use. The work indicator gives a very good indication of how worn the actuator is.</p> <p>Typical minimum lifetime performance without abuse of the actuator</p> <p>LA23: 10.000 cycles in life test equals: 1.500.000 A*S LA31: 10.000 cycles in life test equals: 2.900.000 A*S LA34: 10.000 cycles in life test equals: 4.200.000 A*S LA43: 10.000 cycles in life test equals: 3.750.000 A*S LA44: 10.000 cycles in life test equals: 5.600.000 A*S</p> <p>The indicator can only be reset by exchanging the actuator through Service Data Tool "Advanced Settings" menu. Please contact the lift manufacturer in order to decide when it is appropriate to consider exchanging the actuator.</p>
Cycles  4	<p>Total cycles:</p> <p>The total number of cycles made by the actuator on channel 1.</p> <p>One cycle is defined as; Driving with load (the actuator draws more than 1.5 Amps). Driving direction up for a minimum of 5 seconds (several activations are allowed), followed by driving down for a minimum of 2 seconds.</p> <p>The indicator can only be reset by exchanging the actuator through Service Data Tool "Advanced Settings" menu or by pressing the up and down button on channel 1 for 5 seconds.</p>
Cutoff limit In Out [A] 4,0 4,0	<p>Cut-off limit:</p> <p>To ensure that the lift will stop if the current draw exceeds the preset limits.</p> <p>The best way to set the cut-off limit is via learn mode.</p>



How to conclude on the service information


Total cycles and total work:

Please contact the lift manufacturer in order to decide when it is appropriate to consider exchanging the actuator.

Overload:


If overload has occurred it is recommended to consider stronger lifting equipment with higher working load for the particular patients/institution.

Control box view



Symbol	Data description and recommendations
	<p>Read out how many times a hand control button or the button on the control box has been activated on CH1.</p>
	<p>Battery state:</p> <p>Green; battery fully charged (approx. 100-50% capacity remaining)</p> <p>Yellow; battery approx. half charged (approx. 50-25% capacity remaining)</p> <p>Orange; Battery empty (incl. audio signal when the hand control is activated). Limited actuator function.</p> <p>Good battery behaviour is charging the battery at all times when it is not in use. Deep discharging of the battery can cause damage of the battery.</p>
	<p>Since last reset:</p> <p>Number of times the control box has detected a critical battery level. If batteries are switched between lifts, it does not say anything about the battery condition. It is only an indication of the battery behaviour. The "Battery low critical level" counter will be reset after each service reset, but it will be saved if making a report.</p> <p>Total Number of times the control box has detected a critical battery level. If batteries are switched between lifts, it does not say anything about the battery condition. It is only an indication of the battery behaviour.</p> <p>The number of critical battery levels will not be valid if the emergency stop button is used as an on/off power button and not for emergency situations only.</p>

Controls view



If a diode (or hand control button) is activated, one of the coded buttons will light up as the example above shows.




Codes	Code explanation
<p>CBJC keys</p> <p>1 U 1 D 2 U 2 D</p> <p><input checked="" type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/></p>	<p>Up and down buttons on the CBJC control box. HB40 like signals</p>
<p>Openbus "H-Frame" Input</p> <p>H0 H1 H2 H3 H4 H5 H6 H7 H8 H9 H10 H11 H12 H13 H14 H15 <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/></p> <p>H16 H17 H18 H19 H20 H21 H22 H23 H24 H25 H26 H27 H28 H29 H30 H31 <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/></p>	<p>OpenBus signals</p>
<p>Openbus "A-Frame" Led's</p> <p>A0 A1 A2 A3 A4 A5 A6 A7 A8 A9 A10 A11 <input checked="" type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/></p>	<p>OpenBus Diode signals</p>
<p>Openbus "A-Frame" Input</p> <p>A12 A13 A14 A15 A16 A17 A18 A19 A20 A21 <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/></p> <p>A22 A23 A24 A25 A26 A27 A28 A29 A30 A31 <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/> <input type="radio"/></p>	<p>OpenBus button signals</p>



Report view




After each service visit it is recommended to fill in the service report and press “save” to maintain service data.

Data	Recommended procedure and explanation
<p>Date: <input type="text" value="20. april 2012"/></p> <p>Service period set to: <input type="text" value="360 Days"/> <input type="text" value="8000 Cycles"/></p> <p>Control Box serial number: <input type="text" value="4294967295-4294967295"/></p>	<p>These data are filled in automatically.</p>
<p>Inspected by: <input type="text"/></p> <p>Company: <input type="text"/></p> <p>Lift ID: <input type="text"/></p>	<p>Fill in these data.</p>
<p>Actions taken: (defects, wear or damage) <input type="text"/></p>	<p>Description made in this box will be read out when the report is saved.</p>
<p>Notes: This service visit</p> <p><input type="text" value="003 - Actuator 3 replaced"/></p> <p><input type="text" value="004 - Handset replaced"/></p> <p><input type="text" value="None"/></p> <p><input type="text" value="None"/></p>	<p>Choose up to 4 notes after the service visit. Messages for the next service visit - will be stored and readable at the next visit.</p>
<p>Notes: This service visit</p> <p><input type="text" value="003 - Actuator 3 replaced"/></p> <p><input type="text" value="004 - Handset replaced"/></p>	<p>Up to 4 notes readable from the previous service visit</p>
	<p>Check the lift extra well on the marked spots. Action is recommended in case of noteworthy observations.</p>
	<p>When the lift has been checked and the information completed, press “save” to save as html or a comma separated file. Relevant actuator statistics and system data will be included in the report.</p>
	<p>Resetting the service counters:</p> <ul style="list-style-type: none"> - Days since last service - Cycles since last service <p>Resetting of service is made by pressing 2 buttons on the hand control or the control box at the same time for 5 seconds. The reset buttons are predefined as lifting arm up and lifting arm down. An audio signal will indicate that the timer has been reset.</p>



Advanced settings view

Please note! This section is only intended for trained and authorised service technicians





Data	Recommended procedure and explanation
<p>Control Box Serial number: 4294967295-4294967295</p> <p>Last change of settings: 0 Days ago</p>	<p>Control box serial number: This number is unique and can be found on the label of the control box. It will be updated automatically if the control box is replaced and cannot be changed.</p> <p>Last change of settings: The date of the last programming is shown and cannot be changed, it will be reset when pressing "Save new settings". It will also be reset when service settings are changed via the special LINAK HB.</p>
<p>Service Settings</p> <p>Days between services: 360 days / 12 months</p> <p>Cycles between services: 8000 cycles</p> <p>Value not set list:</p> <ul style="list-style-type: none"> Value not set 180 days / 6 months 240 days / 8 months 300 days / 10 months 360 days / 12 months 420 days / 14 months 480 days / 16 months 540 days / 18 months 600 days / 20 months 660 days / 22 months 720 days / 24 months 	<p>Service settings: The settings are automatically shown, but they can be changed via the drop-down menu.</p> <p>It is also possible to insert days and cycles in service settings which are not available on the drop down list.</p> <p>Only requirement is – the setting of days must be dividable with 30 to calculate the months.</p>
<p>Actuator Statistics</p> <p>Cutoff Limit out: 10.7, 4.0, 4.0</p> <p>Cutoff Limit in: 10.7, 4.0, 4.0</p> <p>Overloads Out: 0, 0, 0</p> <p>Overloads In: 0, 0, 0</p> <p>Work Out: 0, 0, 0</p> <p>Work In: 0, 0, 0</p> <p>Total Cycles: 0</p>	<p>Actuator info: Remember to update actuator info if the control box is exchanged. In this way you maintain the statistical service data information on the actuators.</p> <p>Remember to update "cut-off limit out" and "cut-off limit in" and all other figures if an actuator is replaced.</p>
<p>Get data from CB</p>	<p>Get / retrieve data from the control box:</p> <p>When entering new data, the data will be visible with bold script. "Get data from CB" undo changes and retrieve existing data from the control box.</p>
<p>Save new settings</p>	<p>Remember to save report with existing data before changing settings</p> <p>Save new settings: This will reset the above information with the new data filled in.</p>



Resetting of service interval after service has been carried out

Resetting of service is done by pressing 2 buttons (lifting arm up and lifting arm down on the hand control or on the control box at the same time for 5 seconds. (Buttons that need to be pressed might differ depending on hand control type or control box). After pressing the buttons for 5 seconds, you will receive an audio signal indicating that the timer has been reset. The timer will reset the diodes/clear the display for service symbol, and start counting a new service period.



Info site

The lift manufacturer can add one information site per section and it is accessible via the icon ⓘ

The information site may include further information e.g. service checklist for the service technician to follow, Lift guidance by lift type etc.

Contacts

FACTORIES

Denmark - Headquarters
LINAK A/S
Phone: +45 73 15 15 15
Fax: +45 74 45 80 48
Fax (Sales): +45 73 15 16 13
Web: www.linak.com

China
LINAK (Shenzhen) Actuator Systems, Ltd.
Phone: +86 755 8610 6656
Phone: +86 755 8610 6990
Web: www.linak.cn

Slovakia
LINAK Slovakia s.r.o.
Phone: +421 51 7563 444
Web: www.linak.sk

Thailand
LINAK APAC Ltd.
Phone: +66 33 265 400
Web: www.linak.com

USA
LINAK U.S. Inc.
Americas Headquarters
Phone: +1 502 253 5595
Fax: +1 502 253 5596
Web: www.linak-us.com
www.linak-latinamerica.com

SUBSIDIARIES

Australia
LINAK Australia Pty. Ltd
Phone: +61 3 8796 9777
Fax: +61 3 8796 9778
E-mail: sales@linak.com.au
Web: www.linak.com.au

Austria
LINAK Zweigniederlassung - Österreich (Wien)
Phone: +43 (1) 890 7446
Fax: +43 (1) 890 744615
E-mail: info@linak.de
Web: www.linak.at - www.linak.hu

Belgium
LINAK Actuator-Systems NV/SA
(Belgium & Luxembourg)
Phone: +32 (0)9 230 01 09
E-mail: beinfo@linak.be
Web: www.linak.be - www.fr.linak.be

Brazil
LINAK Do Brasil Comércio De Atuadores Ltda.
Phone: +55 (11) 2832 7070
Fax: +55 (11) 2832 7060
E-mail: info@linak.com.br
Web: www.linak.com.br

Canada
LINAK Canada Inc.
Phone: +1 502 253 5595
Fax: +1 416 255 7720
E-mail: info@linak.ca
Web: www.linak-us.com

Czech Republic
LINAK C&S s.r.o.
Phone: +42 058 174 1814
Fax: +42 058 170 2452
E-mail: info@linak.cz
Web: www.linak.cz - www.linak.sk

Denmark - International
LINAK International
Phone: +45 73 15 15 15
E-mail: info@linak.com
Web: www.linak.com

Denmark - Sales
LINAK Danmark A/S
Phone: +45 86 80 36 11
Fax: +45 86 82 90 51
E-mail: linak@linak-silkeborg.dk
Web: www.linak.dk

Finland
LINAK OY
Phone: +358 10 841 8700
E-mail: linak@linak.fi
Web: www.linak.fi

France
LINAK France E.U.R.L
Phone: +33 (0) 2 41 36 34 34
Fax: +33 (0) 2 41 36 35 00
E-mail: linak@linak.fr
Web: www.linak.fr

Germany
LINAK GmbH
Phone: +49 6043 9655 0
Fax: +49 6043 9655 60
E-mail: info@linak.de
Web: www.linak.de

India
LINAK A/S India Liaison Office
Phone: +91 120 4531797
Fax: +91 120 4786428
E-mail: info@linak.in
Web: www.linak.in

Italy
LINAK ITALIA S.r.l.
Phone: +39 02 48 46 33 66
Fax: +39 02 48 46 82 52
E-mail: info@linak.it
Web: www.linak.it

Japan
LINAK K.K.
Phone: 81-45-533-0802
Fax: 81-45-533-0803
E-mail: linak@linak.jp
Web: www.linak.jp

Malaysia
LINAK Actuators Sdn. Bhd.
Phone: +60 4 210 6500
Fax: +60 4 226 8901
E-mail: info@linak-asia.com
Web: www.linak.my

Netherlands
LINAK Actuator-Systems B.V.
Phone: +31 76 5 42 44 40 /
+31 76 200 11 10
E-mail: info@linak.nl
Web: www.linak.nl

New Zealand
LINAK New Zealand Ltd
Phone: +64 9580 2071
Fax: +64 9580 2072
E-mail: nzsales@linak.com.au
Web: www.linak.com.au

Norway
LINAK Norge AS
Phone: +47 32 82 90 90
E-mail: info@linak.no
Web: www.linak.no

Poland
LINAK Polska
LINAK Danmark A/S (Spółka Akcyjna)
Phone: +48 22 295 09 70 /
+48 22 295 09 71
E-mail: info@linak.pl
Web: www.linak.pl

Republic of Korea
LINAK Korea Ltd.
Phone: +82 2 6231 1515
Fax: +82 2 6231 1516
E-mail: info@linak.kr
Web: www.linak.kr

Slovakia
LINAK Slovakia S.R.O.
Phone: +421 51 7563 444
Web: www.linak.sk

Spain
LINAK Actuadores, S.L.U
Phone: +34 93 588 27 77
Fax: +34 93 588 27 85
E-mail: esma@linak.es
Web: www.linak.es

Sweden
LINAK Scandinavia AB
Phone: +46 8 732 20 00
Fax: +46 8 732 20 50
E-mail: info@linak.se
Web: www.linak.se

Switzerland
LINAK AG
Phone: +41 43 388 31 88
Fax: +41 43 388 31 87
E-mail: info@linak.ch
Web: www.linak.ch - www.fr.linak.ch
www.it.linak.ch

Taiwan
LINAK (Shenzhen) Actuator systems Ltd.
Taiwan Representative office
Phone: +886 2 272 90068
Fax: +886 2 272 90096
E-mail: sales@linak.com.tw
Web: www.linak.com.tw

Turkey
LINAK İth. İhr. San. ve Tic. A.Ş.
Phone: +90 312 4726338
Fax: +90 312 4726635
E-mail: info@linak.com.tr
Web: www.linak.com.tr

United Kingdom & Ireland
LINAK UK Limited
Phone: +44 (0)121 544 2211
Fax: +44 (0)121 544 2552
E-mail: sales@linak.co.uk
Web: www.linak.co.uk

DISTRIBUTORS

Argentina
Novotec Argentina SRL
Phone: 011-4303-8989 / 8900
Fax: 011-4032-0184
E-mail: info@novotecargentina.com
Web: www.novotecargentina.com

Colombia
MEM Ltda
Phone: +[57] (1) 334-7666
Fax: +[57] (1) 282-1684
E-mail: servicioalcliente@memlda.com.co
Web: www.mem.net.co

India
Mechatronics Control Equipments India Pvt Ltd
Phone: +91-44-28558484, 85
E-mail: bala@mechatronicscontrol.com
Web: www.mechatronicscontrol.com

Indonesia
PT. Himalaya Everest Jaya
Phone: +6 221 544 8956 /+6 221 544 8965
Fax: +6 221 619 1925
Fax (Sales): +6 221 619 4658
E-mail: hejplastic-div@centrin.net.id
Web: www.hej.co.id

Israel
NetivTech LTD
Phone: +972 55-2266-535
Fax: +972 2-9900-560
Email: info@NetivTech.com
Web: www.netivtech.com

Singapore
Servo Dynamics Pte Ltd
Phone: +65 6844 0288
Fax: +65 6844 0070
E-mail: servodynamics@servo.com.sg

South Africa
Industrial Specialised Applications CC
Phone: +27 011 466 0346
E-mail: gartht@isagroup.co.za
Web: www.isaza.co.za

United Arab Emirates
Mechatronics
Phone: +971 4 267 4311
Fax: +971 4 267 4312
E-mail: mechtron@emirates.net.ae

Terms of use

LINAK® takes great care in providing accurate and up-to-date information on its products. However, the user is responsible for determining the suitability of LINAK products for a specific application. Due to continual development, LINAK products are subject to frequent modifications and changes. LINAK reserves the rights to conduct modifications, updates, and changes without any prior notice. For the same reason, LINAK cannot guarantee the correctness and actual status of imprinted information on its products.

LINAK uses its best efforts to fulfil orders. However, for the reasons mentioned above, LINAK cannot guarantee availability of any particular product at any given time. LINAK reserves the right to discontinue the sale of any product displayed on its website or listed in its catalogues or in other written material created and produced by LINAK, LINAK subsidiaries, or LINAK affiliates. All sales are subject to the 'Standard Terms of Sale and Delivery for LINAK A/S' available on LINAK websites. LINAK and the LINAK logotype are registered trademarks of LINAK A/S. All rights reserved.